

ZCCHC3 Polyclonal Antibody

Catalog No: #30000



Package Size: #30000-1 50ul #30000-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZCCHC3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ZCCHC3 (NP_149080.2).
Other Names	ZCCHC3; C20orf99; zinc finger CCHC-type containing 3
Accession No.	Swiss-Prot#:Q9NUD5NCBI Gene ID:85364
Uniprot	Q9NUD5
GeneID	85364;
Calculated MW	43kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using ZCCHC3 at 1:3000 dilution.



Immunofluorescence analysis of L929 cells using ZCCHC3 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only