

Kctd19 Polyclonal Antibody

Catalog No: #30004

Package Size: #30004-1 50ul #30004-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Kctd19 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IF
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of mouse Kctd19 (NP_001288102.1).
Other Names	KCTD19; BTB/POZ domain-containing protein KCTD19
Accession No.	Swiss-Prot#:Q562E2NCBI Gene ID:279499
Uniprot	Q562E2
GeneID	279499;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

IF□1:50 - 1:200

Images



Immunofluorescence analysis of mouse testis using Kctd19 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only