

MEX3D Polyclonal Antibody

Catalog No: #30007

Package Size: #30007-1 50ul #30007-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MEX3D Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	A synthetic peptide of human MEX3D.
Other Names	MEX3D; MEX-3D; MEX3; OK/SW-cl.4; RKHD1; RNF193; TINO; RNA-binding protein MEX3D
Accession No.	Swiss-Prot#:Q86XN8NCBI Gene ID:399664
Uniprot	Q86XN8
GeneID	399664;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse kidney, using MEX3D at 1:1000 dilution.

Background

Note: This product is for in vitro research use only