

C10orf62 Polyclonal Antibody

Catalog No: #30008

Package Size: #30008-1 50ul #30008-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C10orf62 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human C10orf62 (NP_001009997.2).
Other Names	C10orf62; bA548K23.1; uncharacterized protein C10orf62
Accession No.	Swiss-Prot#:Q5T681NCBI Gene ID:414157
Uniprot	Q5T681
GeneID	414157;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

IHC□1:50 - 1:200

Images



Immunohistochemistry of paraffin-embedded rat testis using C10orf62 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse testis using C10orf62 at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only