

PRKAA1/PRKAA2 Polyclonal Antibody

Catalog No: #30009



Package Size: #30009-1 50ul #30009-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

|                       |                                                 |
|-----------------------|-------------------------------------------------|
| Product Name          | PRKAA1/PRKAA2 Polyclonal Antibody               |
| Host Species          | Rabbit                                          |
| Clonality             | Polyclonal                                      |
| Isotype               | IgG                                             |
| Purification          | Affinity purification                           |
| Applications          | WB,IF                                           |
| Species Reactivity    | Human,Mouse,Rat                                 |
| Immunogen Description | Recombinant protein of human AMPKa1/AMPKa2.     |
| Other Names           | AMPKa1/AMPKa2                                   |
| Accession No.         | Uniprot:Q13131/P54646GenelD:5562/5563           |
| Uniprot               | Q13131/P54646                                   |
| GenelD                | 5562/5563                                       |
| Calculated MW         | Refer to figures                                |
| SDS-PAGE MW           | 62KDa                                           |
| Formulation           | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage               | Store at -20°C. Avoid freeze / thaw cycles.     |


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100


Images




Western blot analysis of extracts of various cell lines, using AMPKa1/AMPKa2 antibody.



Immunofluorescence analysis of NIH/3T3 cells using AMPKα1/AMPKα2 Rabbit pAb.



Immunofluorescence analysis of U2OS cells using AMPKα1/AMPKα2 Rabbit pAb.



Immunofluorescence analysis of PC-12 cells using AMPKα1/AMPKα2 Rabbit pAb.

Note: This product is for in vitro research use only