

# ZMIZ2 Polyclonal Antibody

Catalog No: #30011



Package Size: #30011-1 50ul #30011-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                          |
|-----------------------|--------------------------------------------------------------------------|
| Product Name          | ZMIZ2 Polyclonal Antibody                                                |
| Host Species          | Rabbit                                                                   |
| Clonality             | Polyclonal                                                               |
| Isotype               | IgG                                                                      |
| Purification          | Affinity purification                                                    |
| Applications          | WB                                                                       |
| Species Reactivity    | Human,Mouse,Rat                                                          |
| Immunogen Description | Recombinant fusion protein of human ZMIZ2 (NP_001287888).                |
| Other Names           | ZMIZ2; NET27; TRAFIP20; ZIMP7; hZIMP7; zinc finger MIZ-type containing 2 |
| Accession No.         | Swiss-Prot#:Q8NF64NCBI Gene ID:83637                                     |
| Uniprot               | Q8NF64                                                                   |
| GeneID                | 83637;                                                                   |
| Calculated MW         | 115kDa                                                                   |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.        |
| Storage               | Store at -20°C                                                           |

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ZMIZ2 at 1:1000 dilution.

## Background

ZMIZ2 and ZMIZ1 (MIM 607159) are members of a PIAS (see MIM 603566)-like family of proteins that interact with nuclear hormone receptors. ZMIZ2 interacts with androgen receptor (AR; MIM 313700) and enhances AR-mediated transcription (Huang et al., 2005 [PubMed 16051670]).

---

Note: This product is for in vitro research use only