

DEPDC6 Polyclonal Antibody

Catalog No: #30012

Package Size: #30012-1 50ul #30012-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DEPDC6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human DEPDC6 (NP_073620.2).
Other Names	DEPTOR; DEP.6; DEPDC6; DEP domain-containing mTOR-interacting protein
Accession No.	Swiss-Prot#:Q8TB45NCBI Gene ID:64798
Uniprot	Q8TB45
GeneID	64798;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:100


Images




Immunohistochemistry of paraffin-embedded human colon using DEPDC6 at dilution of 1:200 (40x lens).




Immunohistochemistry of paraffin-embedded human colon carcinoma using DEPDC6 at dilution of 1:200 (40x lens).




Immunohistochemistry of paraffin-embedded human breast cancer using DEPDC6 at dilution of 1:200 (40x lens).




Immunohistochemistry of paraffin-embedded rat bronchus using DEPDC6 at dilution of 1:200 (40x lens).




Immunohistochemistry of paraffin-embedded human stomach using DEPDC6 at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human lung using DEPDC6 at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human gastric cancer using DEPDC6 at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using DEPDC6 at dilution of 1:200 (40x lens).

Background

Note: This product is for in vitro research use only