

TET3 Polyclonal Antibody

Catalog No: #30014

Package Size: #30014-1 50ul #30014-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TET3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human TET3 (NP_001274420.1).
Other Names	TET3; hCG_40738; methylcytosine dioxygenase TET3
Accession No.	Swiss-Prot#:O43151NCBI Gene ID:200424
Uniprot	O43151
GeneID	200424;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only