

SIDT2 Polyclonal Antibody

Catalog No: #30016



Package Size: #30016-1 50ul #30016-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SIDT2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IF
Species Reactivity	Rat
Immunogen Description	Recombinant protein of human SIDT2.
Other Names	SIDT2; CGI-40; SID1 transmembrane family member 2
Accession No.	Swiss-Prot#:Q8NBJ9NCBI Gene ID:51092
Uniprot	Q8NBJ9
GeneID	51092;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

IF□1:50 - 1:100

Background

Note: This product is for in vitro research use only