

Lamin A / C Polyclonal Antibody

Catalog No: #30017



Package Size: #30017-1 50ul #30017-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Lamin A / C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Lamin A / C
Other Names	LMNA; CDCD1; CDDC; CMD1A; CMT2B1; EMD2; FPL; FPLD; FPLD2; HGPS; IDC; LDP1; LFP; LGMD1B; LMN1; LMNC; LMNL1; MADA; PRO1; lamin
Accession No.	Swiss-Prot#:P02545NCBI Gene ID:4000
Uniprot	P02545
GeneID	4000;
Calculated MW	74kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:200 - 1:1000IHC□1:50 - 1:200IF□1:50 - 1:200


Images




Immunofluorescence analysis of C6 cells using Lamin A / C at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunofluorescence analysis of U-2 OS cells using Lamin A / C at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunohistochemistry of paraffin-embedded human mammary cancer using Lamin A / C at dilution of 1:200 (40x lens).




Immunohistochemistry of paraffin-embedded human uterine cancer using Lamin A / C at dilution of 1:200 (40x lens).




Immunohistochemistry of paraffin-embedded human colon using Lamin A / C at dilution of 1:200 (40x lens).




Immunohistochemistry of paraffin-embedded human stomach using Lamin A / C at dilution of 1:200 (40x lens).




Immunohistochemistry of paraffin-embedded human placenta using Lamin A / C at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded rat brain using Lamin A / C at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded mouse spinal cord using Lamin A / C at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using Lamin A / C at dilution of 1:200 (40x lens).

Background

Note: This product is for in vitro research use only