

RBM47 Polyclonal Antibody

Catalog No: #30021



Package Size: #30021-1 50ul #30021-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RBM47 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human RBM47.
Other Names	RBM47; NET18; RNA-binding protein 47
Accession No.	Swiss-Prot#:A0AV96NCBI Gene ID:54502
Uniprot	A0AV96
GeneID	54502;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only