

CAPRIN2 Polyclonal Antibody

Catalog No: #30026



Package Size: #30026-1 50ul #30026-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CAPRIN2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant protein of human CAPRIN2.
Other Names	EEG1; EEG-1; C1QDC1; RNG140
Accession No.	Swiss-Prot#:Q6IMN6NCBI Gene ID:65981
Uniprot	Q6IMN6
GeneID	65981;
Calculated MW	126kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CAPRIN2 at 1:1000 dilution.

Background

The protein encoded by this gene may regulate the transport of mRNA. It may play a role in the differentiation of erythroblasts. Multiple transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only