

HERC4 Polyclonal Antibody

Catalog No: #30027



Package Size: #30027-1 50ul #30027-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HERC4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human HERC4 (NP_056416.2).
Other Names	HERC4; probable E3 ubiquitin-protein ligase HERC4
Accession No.	Swiss-Prot#:Q5GLZ8NCBI Gene ID:26091
Uniprot	Q5GLZ8
GeneID	26091;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using HERC4 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only