

SLA Polyclonal Antibody

Catalog No: #30028



Package Size: #30028-1 50ul #30028-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLA Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human SLA (NP_006739.2).
Other Names	SLA1; SLAP
Accession No.	Swiss-Prot#:Q13239NCBI Gene ID:6503
Uniprot	Q13239
GeneID	6503;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Background

Note: This product is for in vitro research use only