

SLC7A6 Polyclonal Antibody

Catalog No: #30032



Package Size: #30032-1 50ul #30032-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC7A6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human SLC7A6 (NP_003974.3).
Other Names	LAT3; LAT-2; y+LAT-2
Accession No.	Swiss-Prot#:Q92536NCBI Gene ID:9057
Uniprot	Q92536
GeneID	9057;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only