

Acetyl-FOXO1-K262/K265/K274 pAb

Catalog No: #30033



Package Size: #30033-1 50ul #30033-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Acetyl-FOXO1-K262/K265/K274 pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic acetylated peptide around K262/K265/K274 of human FOXO1.
Other Names	FKH1; FKHR; FOXO1A
Accession No.	Swiss-Prot#:_NCBI Gene ID:2308
GeneID	:2308
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma.

Note: This product is for in vitro research use only