

VPS13A Polyclonal Antibody

Catalog No: #30035

Package Size: #30035-1 50ul #30035-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VPS13A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human VPS13A.
Accession No.	Swiss-Prot#:Q96RL7NCBI Gene ID:23230
Uniprot	Q96RL7
GeneID	23230;
Calculated MW	360kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of U-87MG cells, using VPS13A at 1:1000 dilution.

Background

Note: This product is for in vitro research use only