

VPS13B Polyclonal Antibody

Catalog No: #30036



Package Size: #30036-1 50ul #30036-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	VPS13B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human VPS13B.
Accession No.	Swiss-Prot#:Q7Z7G8NCBI Gene ID:157680
Uniprot	Q7Z7G8
GeneID	157680;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Background

Note: This product is for in vitro research use only