

Smad1/Smad5/Smad9 Polyclonal Antibody

Catalog No: #30038



Package Size: #30038-1 50ul #30038-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Smad1/Smad5/Smad9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human Smad1/Smad5/Smad9.
Accession No.	Swiss-Prot#:_NCBI Gene ID:4086/4090/4093
GeneID	4086;4090;4093;
Calculated MW	70kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using Smad1/Smad5/Smad9 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only