

MAP1S Polyclonal Antibody

Catalog No: #30042

Package Size: #30042-1 50ul #30042-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAP1S Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human MAP1S.
Other Names	MAP1S; BPY2IP1; C19orf5; MAP8; VCY2IP-1; VCY2IP1; microtubule-associated protein 1S
Accession No.	Swiss-Prot#:Q66K74NCBI Gene ID:55201
Uniprot	Q66K74
GeneID	55201;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:1000

Background

Note: This product is for in vitro research use only