

GALK2 Polyclonal Antibody

Catalog No: #30051

Package Size: #30051-1 50ul #30051-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GALK2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human GALK2 (NP_002035.1).
Other Names	GK2
Accession No.	Swiss-Prot#:Q01415NCBI Gene ID:2585
Uniprot	Q01415
GeneID	2585;
Calculated MW	50kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse liver, using GALK2 at 1:1000 dilution.

Background

This gene encodes a highly efficient N-acetylgalactosamine (GalNAc) kinase, which has galactokinase activity when galactose is present at high concentrations. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2014]

Note: This product is for in vitro research use only