

# GK2 Polyclonal Antibody

Catalog No: #30052



Package Size: #30052-1 50ul #30052-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                   |
|-----------------------|-------------------------------------------------------------------|
| Product Name          | GK2 Polyclonal Antibody                                           |
| Host Species          | Rabbit                                                            |
| Clonality             | Polyclonal                                                        |
| Isotype               | IgG                                                               |
| Purification          | Affinity purification                                             |
| Applications          | WB,IF                                                             |
| Species Reactivity    | Human,Mouse,Rat                                                   |
| Immunogen Description | Recombinant fusion protein of human GK2 (NP_149991.2).            |
| Other Names           | GKP2, GKTA                                                        |
| Accession No.         | Swiss-Prot#:Q14410NCBI Gene ID:2712                               |
| Uniprot               | Q14410                                                            |
| GeneID                | 2712;                                                             |
| Calculated MW         | 56kDa                                                             |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage               | Store at -20°C                                                    |


## Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100


## Images



Western blot analysis of extracts of various cell lines, using GK2 at 1:1000 dilution.



Immunofluorescence analysis of C6 cells using GK2 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of L929 cells using GK2 at dilution of 1:100. Blue: DAPI for nuclear staining.

## Background

Note: This product is for in vitro research use only