

EXOSC10 Polyclonal Antibody

Catalog No: #30062



Package Size: #30062-1 50ul #30062-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	EXOSC10 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human EXOSC10 (NP_001001998.1).
Other Names	PM-Scl, PM/Scl-100, PMSCL, PMSCL2, RRP6, Rrp6p, p2, p3, p4
Accession No.	Swiss-Prot#:Q01780NCBI Gene ID:5394
Uniprot	Q01780
GeneID	5394;
Calculated MW	110kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using EXOSC10 at 1:1000 dilution.



Western blot analysis of extracts of Rat brain, using EXOSC10 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only