

PNLIPRP1 Polyclonal Antibody

Catalog No: #30063

Package Size: #30063-1 50ul #30063-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PNLIPRP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PNLIPRP1 (NP_006220.1).
Other Names	PLRP1
Accession No.	Swiss-Prot#:P54315NCBI Gene ID:5407
Uniprot	P54315
GeneID	5407;
Calculated MW	52kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:500 - 1:2000IHC□1:100 - 1:200IF□1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded rat pancreas using PNLIPRP1 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human appendix using PNLIPRP1 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded mouse pancreas using PNLIPRP1 at dilution of 1:100 (40x lens).




Western blot analysis of extracts of various cell lines, using PNLIPRP1 at 1:1000 dilution.



Immunofluorescence analysis of rat pancreas using PNLIPRP1 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of rat pancreas using PNLIPRP1 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of mouse pancreas using PNLIPRP1 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only