

SH3BGRL Polyclonal Antibody

Catalog No: #30066



Package Size: #30066-1 50ul #30066-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SH3BGRL Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SH3BGRL (NP_003013.1).
Other Names	HEL-S-115, SH3BGR
Accession No.	Swiss-Prot#:O75368NCBI Gene ID:6451
Uniprot	O75368
GeneID	6451;
Calculated MW	16kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:100 - 1:200IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using SH3BGRL at 1:1000 dilution.



Immunofluorescence analysis of U-2 OS cells using SH3BGRL at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only