

SPRR1A Polyclonal Antibody

Catalog No: #30069



Package Size: #30069-1 50ul #30069-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SPRR1A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of Human SPRR1A.
Other Names	SPRK
Accession No.	Swiss-Prot#:P35321NCBI Gene ID:6698
Uniprot	P35321
GeneID	6698;
Calculated MW	19kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using SPRR1A at 1:1000 dilution.



Immunofluorescence analysis of rat skin using SPRR1A at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of human skin using SPRR1A at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of mouse skin using SPRR1A at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only