

TMEM11 Polyclonal Antibody

Catalog No: #30081



Package Size: #30081-1 50ul #30081-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TMEM11 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TMEM11 (NP_003867.1).
Other Names	C17orf35, PM1, PMI
Accession No.	Swiss-Prot#:P17152NCBI Gene ID:8834
Uniprot	P17152
GeneID	8834;
Calculated MW	22kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images



Western blot analysis of extracts of HeLa cells, using TMEM11 at 1:1000 dilution.



Immunofluorescence analysis of C6 cells using TMEM11 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of L929 cells using TMEM11 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only