

LDB1 Polyclonal Antibody

Catalog No: #30082



Package Size: #30082-1 50ul #30082-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LDB1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human LDB1 (NP_003884.1).
Other Names	CLIM-2, CLIM2, LDB-1, NLI
Accession No.	Swiss-Prot#:Q86U70NCBI Gene ID:8861
Uniprot	Q86U70
GeneID	8861;
Calculated MW	46kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB1:500 - 1:2000IHC1:50 - 1:100


Images




Immunohistochemistry of paraffin-embedded human uterine cancer using LDB1 at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human tonsil using LDB1 at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human colon using LDB1 at dilution of 1:200 (40x lens).



Western blot analysis of extracts of various cell lines, using LDB1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only