

LIMD1 Polyclonal Antibody

Catalog No: #30085



Package Size: #30085-1 50ul #30085-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LIMD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human LIMD1 (NP_055055.1).
Accession No.	Swiss-Prot#:Q9UGP4NCBI Gene ID:8994
Uniprot	Q9UGP4
GeneID	8994;
Calculated MW	72kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using LIMD1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only