

ITGBL1 Polyclonal Antibody

Catalog No: #30087



Package Size: #30087-1 50ul #30087-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ITGBL1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human ITGBL1 (NP_004782.1).
Other Names	OSCP, TIED
Accession No.	Swiss-Prot#:O95965NCBI Gene ID:9358
Uniprot	O95965
GeneID	9358;
Calculated MW	54kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ITGBL1 at 1:1000 dilution.

Background

This gene encodes a beta integrin-related protein that is a member of the EGF-like protein family. The encoded protein contains integrin-like cysteine-rich repeats. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012]

Note: This product is for in vitro research use only