

SUPT7L Polyclonal Antibody

Catalog No: #30090

Package Size: #30090-1 50ul #30090-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SUPT7L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SUPT7L (NP_001269658).
Other Names	SPT7L, STAF65, STAF65(gamma), STAF65G, SUPT7H
Accession No.	Swiss-Prot#:O94864NCBI Gene ID:9913
Uniprot	O94864
GeneID	9913;
Calculated MW	60kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of Rat lung, using SUPT7L at 1:1000 dilution.

Background

SUPT7L is a protein subunit of the human STAGA complex (SPT3; (MIM 602947)/TAF9 (MIM 600822)/GCN5 (MIM 602301) acetyltransferase complex), which is a chromatin-modifying multiprotein complex (Martinez et al., 2001 [PubMed 11564863]).[supplied by OMIM, Apr 2009]

Note: This product is for in vitro research use only