

TOM1L1 Polyclonal Antibody

Catalog No: #30091

Package Size: #30091-1 50ul #30091-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TOM1L1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of Human TOM1L1.
Other Names	SRCASM
Accession No.	Swiss-Prot#:O75674NCBI Gene ID:10040
Uniprot	O75674
GeneID	10040;
Calculated MW	53kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TOM1L1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only