

LPCAT3 Polyclonal Antibody

Catalog No: #30093



Package Size: #30093-1 50ul #30093-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LPCAT3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms,Rt,Hu
Immunogen Description	Recombinant fusion protein of human LPCAT3 (NP_005759.4).
Other Names	C3F, LPCAT, LPLAT 5, LPSAT, MBOAT5, OACT5, nassy
Accession No.	Swiss-Prot#:Q6P1A2NCBI Gene ID:10162
Uniprot	Q6P1A2
GeneID	10162;
Calculated MW	56kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using LPCAT3 at 1:1000 dilution.

Note: This product is for in vitro research use only