

RBM7 Polyclonal Antibody

Catalog No: #30094



Package Size: #30094-1 50ul #30094-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RBM7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human RBM7 (NP_057174.1).
Accession No.	Swiss-Prot#:Q9Y580NCBI Gene ID:10179
Uniprot	Q9Y580
GeneID	10179;
Calculated MW	36kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RBM7 at 1:1000 dilution.



Western blot analysis of extracts of Mouse thymus, using RBM7 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only