

WBP4 Polyclonal Antibody

Catalog No: #30098

Package Size: #30098-1 50ul #30098-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	WBP4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Rat
Immunogen Description	Recombinant fusion protein of human WBP4 (NP_009118.1).
Other Names	FBP21
Accession No.	Swiss-Prot#:O75554NCBI Gene ID:11193
Uniprot	O75554
GeneID	11193;
Calculated MW	43kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Images



Western blot analysis of extracts of Rat thymus, using WBP4 at 1:1000 dilution.

Background

This gene encodes WW domain-containing binding protein 4. The WW domain represents a small and compact globular structure that interacts with proline-rich ligands. This encoded protein is a general spliceosomal protein that may play a role in cross-intron bridging of U1 and U2 snRNPs in the spliceosomal complex A. [provided by RefSeq, Jul 2008]

Note: This product is for in vitro research use only