

MON1B Polyclonal Antibody

Catalog No: #30102



Package Size: #30102-1 50ul #30102-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MON1B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human MON1B (NP_055755.1).
Other Names	HSRG1, SAND2, SRG1
Accession No.	Swiss-Prot#:Q7L1V2NCBI Gene ID:22879
Uniprot	Q7L1V2
GeneID	22879;
Calculated MW	59kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of U-87MG cells, using MON1B at 1:1000 dilution.

Background

Note: This product is for in vitro research use only