

SCG5 Antibody

Catalog No: #46085

Package Size: #46085-1 50ul #46085-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SCG5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SCG5 Antibody detects endogenous levels of total SCG5
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SCG5
Target Name	SCG5
Other Names	SCG5, 7B2, Neuroendocrine protein 7B2, p7B2, Secretogranin V (7B2 protein), Secretogranin-5, SGNE1, SgV, Secretogranin V, Pituitary polypeptide
Accession No.	Swiss-Prot#: P05408NCBI Gene ID:6447
Uniprot	P05408
GeneID	6447;
Calculated MW	24 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only