

CEP131 Polyclonal Antibody

Catalog No: #30105

Package Size: #30105-1 50ul #30105-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CEP131 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of Human CEP131.
Other Names	AZ1, AZI1, ZA1
Accession No.	Swiss-Prot#:Q9UPN4NCBI Gene ID:22994
Uniprot	Q9UPN4
GeneID	22994;
Calculated MW	145kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of HeLa cells, using CEP131 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only