

LIMCH1 Polyclonal Antibody

Catalog No: #30106

Package Size: #30106-1 50ul #30106-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LIMCH1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human LIMCH1 (NP_001106189.1).
Other Names	LIMCH1A, LMO7B
Accession No.	Swiss-Prot#:Q9UPQ0NCBI Gene ID:22998
Uniprot	Q9UPQ0
GeneID	22998;
Calculated MW	140kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of LO2 cells, using LIMCH1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only