

ANKLE2 Polyclonal Antibody

Catalog No: #30108



Package Size: #30108-1 50ul #30108-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ANKLE2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ANKLE2 (NP_055929.1).
Other Names	KIAA0692;LEMD7;Lem4;MCPH16;ANKLE2
Accession No.	Uniprot:Q86XL3GeneID:23141
Uniprot	Q86XL3
GeneID	23141
Calculated MW	104kDa
SDS-PAGE MW	104kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using ANKLE2 antibody.



Immunohistochemistry of paraffin-embedded Mouse kidney using ANKLE2 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Human lung cancer using ANKLE2 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Rat testis using ANKLE2 Rabbit pAb.

Background

This gene encodes a member of the LEM family of inner nuclear membrane proteins. The encoded protein functions as a mitotic regulator through postmitotic formation of the nuclear envelope. Mutations in this gene cause morphology defects in the nuclear envelope and BAF hyperphosphorylation. [provided by RefSeq, Mar 2014]

Note: This product is for in vitro research use only