

FAM120A Polyclonal Antibody

Catalog No: #30113



Package Size: #30113-1 50ul #30113-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FAM120A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human FAM120A (NP_055427.2).
Other Names	C9orf10, HBVPTPAP, OSSA
Accession No.	Swiss-Prot#:Q9NZB2NCBI Gene ID:23196
Uniprot	Q9NZB2
GeneID	23196;
Calculated MW	122kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using FAM120A at 1:1000 dilution.

Background

Note: This product is for in vitro research use only