

ZNF500 Polyclonal Antibody

Catalog No: #30118

Package Size: #30118-1 50ul #30118-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF500 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human ZNF500 (NP_067678.1).
Other Names	ZKSCAN18, ZSCAN50
Accession No.	Swiss-Prot#:O60304NCBI Gene ID:26048
Uniprot	O60304
GeneID	26048;
Calculated MW	54kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse brain, using ZNF500 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only