

FAM127B Polyclonal Antibody

Catalog No: #30119



Package Size: #30119-1 50ul #30119-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FAM127B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human FAM127B (NP_001071640.1).
Other Names	CXX1b, FAM127B, MAR8A, SIRH6
Accession No.	Swiss-Prot#:Q9BWD3NCBI Gene ID:26071
Uniprot	Q9BWD3
GeneID	26071;
Calculated MW	18kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using FAM127B at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using FAM127B at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only