

GMIP Polyclonal Antibody

Catalog No: #30134



Package Size: #30134-1 50ul #30134-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GMIP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human GMIP (NP_057657.2).
Other Names	ARHGAP46;GMIP
Accession No.	Uniprot:Q9P107GeneID:51291
Uniprot	Q9P107
GeneID	51291
Calculated MW	107kDa
SDS-PAGE MW	107kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using GMIP antibody.



Immunohistochemistry of paraffin-embedded Mouse spleen using GMIP Rabbit pAb.



Immunohistochemistry of paraffin-embedded Human liver cancer using GMIP Rabbit pAb.



Immunohistochemistry of paraffin-embedded Rat spleen using GMIP Rabbit pAb.

Background

This gene encodes a member of the ARHGAP family of Rho/Rac/Cdc42-like GTPase activating proteins. The encoded protein interacts with the Ras-related protein Gem through its N-terminal domain. Separately, it interacts with RhoA through a RhoGAP domain, and stimulates RhoA-dependent GTPase activity. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]

Note: This product is for in vitro research use only