

SETD4 Polyclonal Antibody

Catalog No: #30139



Package Size: #30139-1 50ul #30139-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SETD4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human SETD4 (NP_001007260).
Other Names	C21orf18, C21orf27
Accession No.	Swiss-Prot#:Q9NVD3NCBI Gene ID:54093
Uniprot	Q9NVD3
GeneID	54093;
Calculated MW	50kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using SETD4 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only