

NDC1 Polyclonal Antibody

Catalog No: #30149



Package Size: #30149-1 50ul #30149-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NDC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NDC1 (NP_060557.3).
Other Names	NET3, TMEM48
Accession No.	Swiss-Prot#:Q9BTX1NCBI Gene ID:55706
Uniprot	Q9BTX1
GeneID	55706;
Calculated MW	76kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NDC1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only