

UTP6 Polyclonal Antibody

Catalog No: #30151



Package Size: #30151-1 50ul #30151-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	UTP6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human UTP6 (NP_060898.2).
Other Names	C17orf40, HCA66
Accession No.	Swiss-Prot#:Q9NYH9NCBI Gene ID:55813
Uniprot	Q9NYH9
GeneID	55813;
Calculated MW	70kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using UTP6 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only