

PNO1 Polyclonal Antibody

Catalog No: #30153



Package Size: #30153-1 50ul #30153-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PNO1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human PNO1 (NP_064528.1).
Other Names	KHRBP1, RRP20
Accession No.	Swiss-Prot#:Q9NRX1NCBI Gene ID:56902
Uniprot	Q9NRX1
GeneID	56902;
Calculated MW	28kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PNO1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only