

HHATL Polyclonal Antibody

Catalog No: #30155

Package Size: #30155-1 50ul #30155-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HHATL Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human HHATL (NP_065758.3).
Other Names	C3orf3, GUP1, MBOAT3, MSTP002, OACT3
Accession No.	Swiss-Prot#:Q9HCP6NCBI Gene ID:57467
Uniprot	Q9HCP6
GeneID	57467;
Calculated MW	57kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Rat brain, using HHATL at 1:1000 dilution.

Background

Note: This product is for in vitro research use only