

EEFSEC Polyclonal Antibody

Catalog No: #30160



Package Size: #30160-1 50ul #30160-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	EEFSEC Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human EEFSEC (NP_068756.2).
Other Names	EFSEC, SELB
Accession No.	Swiss-Prot#:P57772NCBI Gene ID:60678
Uniprot	P57772
GeneID	60678;
Calculated MW	72kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using EEFSEC at 1:1000 dilution.

Background

Note: This product is for in vitro research use only